GGT BR

IFCC


Enzymatic colorimetric method

KINETIC

REF CT10192 3 x 50 mL CONTENTS R1.Reagent 3 x 40 mL R2.Reagent 3 x 10 mL	REF CT10195 12 x 50 mL CONTENTS R1.Reagent 12 x 40 mL R2.Reagent 12 x 10 mL
For <i>in vitro</i> diagnostic use only	

SUMMARY

Gamma-glutamyltransferase (γ -GT) catalyzes the transfer of a γ -glutamyl group from γ -glutamyl-3-carboxy-4-nitroanilide to glycylglycine with the formation of L- γ -glutamyl-glycylglycine and 5 amino-2-nitro-benzoate. The amount of 5-amino-2-nitro-benzoate formed, monitored kinetically at 405 nm, is proportional to the enzyme activity present in the sample.¹

REAGENTS

R1. Buffer. TRIS 133 mmol/L, glycylglycine 138 mmol/L.

R2. Substrate. L- γ -Glutamyl-3-carboxy-4-nitroanilide 23 mmol/L.

PREPARATION

The Reagents are ready-to-use.

STORAGE AND STABILITY

Store at 2-8°C. The Reagents are stable until the expiry date stated on the label. After daily use stored tightly closed and protected from light. On board the reagents are stable 30 days. Discard the reagent if the blank presents an absorbance over 1.200 at 405 nm against distilled water.

SAMPLE COLLECTION

Serum or EDTA plasma free of hemolysis. Fluoride, citrate and oxalate inhibit γ -GT activity². The enzyme in the sample is stable for at least 1 week at 2-8°C and for at least 2 months when frozen.

INTERFERENCES⁶

- Lipemia (intralipid >2.5 g/L) may affect the results.
- Bilirubin (> 10 mg/dL) may affect the results.
- Hemoglobin (> 8 g/L) may affect the results.
- Other drugs and substances may interfere⁶.

INSTRUMENTATION AND MATERIALS

- LIDA analyzer.
- Laboratory equipment.
- iso-Clean Solution. Ref. CT18002.
- Multicalibrator CC/H 10x3 mL Ref. CT19750

AUTOMATED PROCEDURE

A graphic display pictures the specific sets corresponding to the technical application outlined for this test. Any new application, to the instrument should be validated to confirm that results meet the analytical performance of the method. It is recommended to validate periodically the instrument.

CALIBRATION

No calibration is required. Activity is calculated against theoretical factor.

CALCULATION

Patient values are calculated automatically by theoretical factor.

RESULTS

Samples with $\Delta A/\text{min}$ exceeding 0.200 at 405 nm should be diluted 1:10 with saline and assayed again. Multiply the results by 10.

If results are to be expressed as SI units apply: $U/L \times 16.67 = \text{nkat/L}$

EXPECTED VALUES⁴

Serum, plasma

Temperature	37°C	30°C	25°C
Men	10-50 U/L (167-834 nKat/L)	7-35 U/L (117-538 nKat/L)	5-25 U/L (83-417 nKat/L)
Women	8-35 U/L (133-583 nKat/L)	6-25 U/L (100-417 nKat/L)	6-25 U/L (83-300 nKat/L)

It is recommended that each laboratory establishes its own reference range.

QUALITY CONTROL

To ensure adequate quality control (QC), each run should include a set of controls (normal and abnormal) with assayed values handled as unknowns. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

DIAGNOSTIC CHARACTERISTICS

Gamma-glutamyl transferase is the most sensitive enzymatic indicator of hepatobiliary disease available. It is highest in cases of *intrahepatic* or *posthepatic biliary obstruction*, and is more sensitive than alkaline phosphatase in detecting *obstructive jaundice*, *colangitis*, and *cholecystitis*.

Elevated levels are noted in the sera of patients with *alcoholic cirrhosis* and from people who are heavy drinkers. The enzyme levels are important in detecting alcohol induced liver disease correlating well with the duration of the drug action.

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Performance characteristics are available on request.

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In vitro diagnostic reagent for the determination of gamma-glutamyltransferase (γ -GT) in serum or plasma.

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Factor included in the application a 37°C

CONTROL CALIDAD *QUALITY CONTROL*

Linear Human Multisera Normal CT19800

Linear Human Multisera Abnormal CT19850

DETERMINACIONES *NUMBER OF TESTS*

600 tests/kit, Code CT10192

2.400 tests/kit, Code CT10195

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(dead volume is not taken in consideration)

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For additional information read reagent packaging insert.

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The reported application is given only as a guideline. Any application to an instrument should be validated to demonstrate that results meet the performance characteristics of the method.

CALIBRATION PARAMETERS

Rule:	One-point linear (K Factor)	Calibrator 1:	
Factor K	1158	Calibrator 2:	
Replicates	1	Calibrator 3:	
Interval		Calibrator 4:	
Sensitivity		Calibrator 5:	
Coefficient		Calibrator 6:	
Error Limit			
Blank response			
Difference Limit			
SD			

TEST PARAMETERS

No.		
Test	GGT	
Method	Kinetic	
Direction	Ascend	
Unit	U/L	
Decimals	2	
Prim. Wave.	405	
Sec. Wave.	0	
Sample Vol. (µL)	25	
R1 Vol. (µL)	200	
R2 Vol. (µL)	50	
Line. Limit	20	
Incubation	10	
Reaction	6	19
Antigen	☐ Check Antigen	
Substrate	0	
Response	L Limit 0	U Limit 0
S. Volume	0	
Ratio	0	
R1 Blank	L Limit 0	U Limit 0
Mix. R Blank	L Limit 0	U Limit 0
Linearity	L Limit 0	U Limit 800
Full Name	GGT	
Print No.		

For other instrument software version please contact with LINEAR.